



CX1. Spinerogy® wheel with 18 black spokes, complete with tyre and pushrim

Spinerogy® wheels

Different wheel combinations

Spinerogy wheel are made of double-wall rim in aluminium with 12, 18 or 24 fibre spokes mounted in a 1.87" or 2.3" aluminium hub. Extremely reliable and comfortable for your everyday needs. Very quiet rolling wheels (no spoke creaking). Ultra lightweight and attractive design. Minimal maintenance required.

We have on stock the Spinerogy® wheels Everyday, LX or SLX in different sizes and colours (see table below). As standard we supply the Spinerogy® wheels with black spokes. On request other spoke colours are also available: black, yellow, red, blue or white (longer delivery times possible).

Available in sizes of 24 ", 25" and 26"

Spinerogy Hub types



CX2a. Spinerogy standard hub 1,80", black or silver coloured



CX2b. Spinerogy standard hub, 1,80" CNC machined, black or silver



CX2c. Spinerogy large hub, 2,30", black or silver coloured

	hub width	hub type	image	Ø of spokes	no. of spokes
Everyday	1,80" - 46,5 mm	standard	CX2a	3 mm	18
LX	1,80" - 46,5 mm	CNC machined	CX2b	4 mm	12
SLX	2,30" - 60 mm	large	CX2c	4 mm	24



Spinerogy with yellow spokes



Spinerogy with red spokes



Spinerogy with blue spokes



Spinerogy with white spokes

Tires and push-rims are available as an option!

On request the Spinerogy® wheels can be supplied with mounted tires and pushrims. Various types of tires and push-rims are available (see the price list for details).



CX4. "Ultralite" wheel

The **Ultralite** wheel is an easy and inexpensive wheel with a total weight of only **1120 grams** including tire, tube and grip ring. The load capacity is 75 kg.

The standard version comes with double walled black rim, black aluminium hub, aluminium push-rim and 1" tires.

The wheel can be delivered with 12 mm or with 1/2" shafts. It is only available on 24".



CX5. "Off-Road" wheel with mountain-bike tyre

24" **Off-road** wheel, complete with black double wall rim, asymmetrically mounted on black high flange hub, with black aluminium push-rim, with tube and black 24" x 1.75" off-road tyre.

Available with 12 mm or 1/2" shafts.