



CX1. Spinerogy® SPOX Everyday wheel with 18 black spokes Complete with tyre/tube and pushrim.

Spinerogy® wheel

Different wheel combinations

We have in stock Spinerogy wheels Everyday, Sport, LX or SLX in different sizes and colours

Spinerogy wheel are made of double-wall rim in aluminium with 12, 18 or 24 fiber spokes mounted in a 1.87" or 2.3" aluminium hub. The modern design is extreme robust for influence and give a noiseless drive because of the fibre spokes.

On stock we have Spinerogy wheels with the following specifications:

	Hub distance	Hub type	Ø Spokes	Spokes quantity
SPOX Everyday	1,87" - 47 mm	Standard	3 mm	18
SPOX Sport	2,30" - 58 mm	Standard	3 mm	24
Sport LX	1,87" - 47 mm	Machined	4 mm	12
Sport SLX	2,30" - 58 mm	Machined	4 mm	18



SPOX with yellow spokes



SPOX with red spokes



SPOX with blue spokes



SPOX with white spokes

Spinerogy wheels complete with tyre and pushrim

We can deliver against more charge Spinerogy wheels with tyre/tube and pushrim. See pricelist for tyre and pushrim to see the different possibilities.



CX4. "Octopus" fibre wheel

The Octopus fibre wheel is a light weight wheel with a weight inclusive tyre/tube and push-rim of only 1420 grammes. The max. load is 75 kg. It is standard supplied with silver polished aluminium hub and quick release axle with push button, double walled black rim, black fibre spokes, black aluminium push-rim and 1" tyre.

The fibre spokes guarantees strength and a long life-cycle

The wheel can be supplied as well with 12 mm hub as 1/2" hub. We stock 24" and 25", but the wheel is also available in 26" and 28"



CX5. "Tristar" wheel

The Tristar wheel is an aluminium wheel with 3 spokes. The weight incl. tyre/tube and push-rim is 1730 grammes. The colour is black with dots called "black pearl". The max. load is 75 kg. It is standard supplied with quick release axle with push button, double walled black rim, black fibre spokes, black aluminium push-rim and 1" tyre. The strong welded construction guarantees strength and a long life-cycle

The wheel can be supplied as well with 12 mm hub as 1/2" hub. We stock 24" but the wheel is also available in 22".



CX6. "X-Core" wheel

X-Core wheel is made from carbon fibre, and has a max. load of 115 kg. The wheel incl. tyre and push-rim has a weight 1800 grammes.

Is standard supplied with 1" tyre and black aluminium push-rim

X-Core is only available in 25". Can be supplied as well with 12 mm as 1/2" hub.



CX4 Typhoon wheel

The Typhoon fibre wheel has a wide hub with a more "raw" look. The weight is 1470 grammes incl. tyre/tube and push-rim The max. load is 75 kg. It is standard supplied with silver/black aluminium hub and quick release axle with push button, double walled black rim, black fibre spokes, black aluminium push-rim and 1" tyre. The fibre spokes guarantees strength and a long life-cycle


The wheel can be supplied as well with 12 mm hub as 1/2" hub. We stock 24".



CX8. "Ultralite" wheel

The Typhoon fibre wheel has a wide hub with a more "raw" look. The weight is 1470 grammes incl. tyre/tube and push-rim The max. load is 75 kg. It is standard supplied with silver/black aluminium hub and quick release axle with push button, double walled black rim, black fibre spokes, black aluminium push-rim and 1" tyre. The fibre spokes guarantees strength and a long life-cycle

The wheel can be supplied as well with 12 mm hub as 1/2" hub. We stock 24".



CX9. Off-Road wheel

24" Off-road wheel, complete with black double wall rim, asymmetrically mounted on black high flange hub, with black aluminium push-rim, with tube and black 24" x 1.75" off-road tyre.

Available with 12 mm or 1/2" shafts.