

### Detachable toilet seats



**E1.** Raised toilet seat without lid



**E2.** Raised toilet seat with lid

The ergonomic raised toilet seat is available in heights of 5, 10 or 15 cm.

The toilet seat opening is 220 mm wide and 270 mm deep. The maximum load is 190 kg.

Since no edges are present, where dirt can accumulate, the raised toilet seat is easy to clean. It is waterproof and can be sterilized at 85°C.

### Fixed raised toilet seats



**E3a.** Fixed raised toilet seat

#### Fixed 10 cm toilet seat

Fixed toilet seat (E3) with lid. New design with new developed fitting system in aluminium. The Toilet seat is easy to fit in the holes where the original seat was fitted.



**E3b.** Fitting system

The armrests can be lifted up and are also detachable.

The opening is 220 mm wide and 270 mm deep.

Max. load is 120 kg

#### Height adjustable fixed toilet seat

The height adjustable toilet seat (E4) is made from non-corroding materials and is easy to clean. Can be fitted on all standards toilets. Has lift-up and detachable armrests that can be taken off in either side.

Max. load 120 kg

Adjustable in heights to 80 mm, 120 mm and 160 mm  
Width inclusive armrests: 610 mm  
Toilet opening is 230 mm x 280 mm (W x D)



**E4a.** Height adjustable toilet seat

### Soft seats



**E5a/E5b.** Soft seats for standard toilet seats (black, white)

The polyurethane-seats (E5) fits standard toilet seats, and can be cleaned with hot water and neutral soap. Available in black and white.

Dimensions: 37.5 x 43.5 cm (W x D)  
The toilet opening is 22 x 27.5 cm (W x D)  
Thickness 4 cm

### Spare-parts for raised toilet seats



**E1c.** Locking device for raised toilet seats 3011..



**E3c.** Fixation kit for fixed raised toilet seat 301701



**E4b.** Fixation kit for fixed raised toilet seat 301702