

Raised toilet seats



E1. Raised toilet seat



E2. Raised toilet seat with lid

An ergonomic designed toilet seat available in heights of 5,10 and 15cm. The opening is 220 mm wide and 270 mm deep. Max load is 190kg. It is easy to clean, as there are no edges where dirt can stick. Is waterproof and can be sterilized at 85°C.

Fixed raised toilet seat.



E3a. Fixed raised toilet seat

Fixed 10 cm toilet seat

Fixed toilet seat (E3) with lid. New design with new developed fitting system in aluminium. The Toilet seat is easy to fit in the holes where the original seat was fitted. The armrests can be lifted up and are also detachable.



E3b. Fitting system

The opening is 220 mm wide and 270 mm deep.

Max. load is 120 kg



E4. Height adjustable toilet seat

Height adjustable fixed toilet seat

The height adjustable toilet seat (E4) is made from non-corroding materials and is easy to clean. Can be fitted on all standards toilets. Has lift-up and detachable armrests that can be taken off in either side.

Max. load 120 kg

Adjustable in heights to: 80 mm, 120 mm and 160 mm

Width inclusive armrests: 610 mm

Toilet opening is (W x D)...: 230 mm x 280 mm



E5. Soft seat for standard toilet seats

Soft seat

The polyurethane-seats (E5) fits standard toilet seats, and can be cleaned with hot water and neutral soap. Available in black and white.

Dimension: 37.5 x 43.5 cm (W x D)

Opening ..: 22 x 27.5 cm (W x D)

Thickness : 4 cm