



**G1.** Anti-decubitus mattress

### Standard Model: Type 418 and 820

Tubular mattress from PU/nylon with double chamber system with alternating pressure. Good standard mattresses for an affordable price. Standard supplied with a mechanical pump, that are low-noise and can be adjusted according to the weight of the patient. The pump is quiet and easy to adjust. A warning light gives warning in case of low air pressure due to a leak in the mattress.

Cycles time is approx. 10 min.

#### Supplied in the dimensions:

Model 418: 80 x 190 x 10 cm (W x L x H)

Model 820: 80 x 190 x 20 cm (W x L x H)

### Comfort Model: Type 516

Tubular mattress from Nylon/PU with double chamber system with alternating pressure. A very good material for high comfort as it is soft and elastic. Supplied with a mechanical pump, that are low-noise and can be adjusted according to the weight of the patient. The pump is quiet and easy to adjust. A warning light gives warning in case of low air pressure due to a leak in the mattress.

Cycles time is approx. 10 min.

#### Supplied in the dimension:

90 cm x 200 cm x 13 cm (W x L x H)

For both mattresses is included an elastic, waterproof cover, made from anti-bacterial material. It has a friction surface to prevent the patient from sliding.



**G2.** Pressure relieving mattress

### Low-cost model: Type 250E

Mattress from the bobble type with double chamber system.

The mattress is made from anti-bacterial, waterproof material with friction cover preventing the patient from sliding.

Included are a mechanical pump where the air pressure can be regulated according to the weight of the patient. The pump is quiet and easy to adjust. A warning light gives warning in case of low air pressure due to a leak in the mattress.

Cycles time is approx. 6 min.

#### Supplied in the dimension:

90 cm x 200 cm x 6 cm (W x L x H)