

## L1. Channel-ramps - CH-ramps

The ramp is made of aluminium with a friction surface. The outside dimensions are 23.5 x 180 cm. Consists of two separate drive-on rails that are hinged at the middle. Folded the ramp can be carried with 2 handles like a suitcase. Used as pairs. Max. load 300 kg.



**L1a.** CH Ramps in use



**L1c.** The CH ramps are hinged in the middle



**L1b.** CH Ramps folded for transport

## L2. Telescopic ramps type TR and LR

The ramps are made of aluminium and with a friction surface. They consist of two separate drive-on rails with 3 segments, which are provided with a locking mechanism to prevent inadvertent extension and retraction of the segments. The ramps are equipped with handles. Used as pairs.

The telescopic ramps are available in different sizes:

- Type TR: 20 x 150 cm, 20 x 210 cm and 20 x 300 cm
- Type LR: 11 x 150 cm and 11 x 210 cm

Max. load 225 kg for the TR type, 180 kg for the LR type.

The LR types are equipped with a handy carrying bag that can be attached to the wheelchair. The carrying bag is optionally also available for the models TR types as an accessory (extra charge).



**L2a.** Telescopic ramps in use



**L2c.** Detail of the locking mechanism



**L2b.** Telescopic ramps telescoped for transport

### L3. Single-folded ramps - SF-ramps

The ramps are made of aluminium with a friction surface and are collapsible in the longitudinal direction - a so-called flat-ramp. This offers the opportunity for volunteers to be able to stand on the ramp together with the wheelchair and the patient. The ramps are equipped with a carrying handle and hinged at the middle they are portable like a suitcase.

Available in 4 sizes : 71 x 40 cm, 76 x 90 cm, 76 x 120 cm and 76 x 150 cm.  
Max load 300 kg .



L3a. SF-040 ramp in use



L3b. SF-040 ramp folded for transport



L3c. SF-ramp in use



L3d. SF ramp folded for transport

### L4. Single-folded ramps - SC-ramps

Inexpensive suitcase ramp in aluminium with anti slip surface.  
Hinged at the middle and carried like a suitcase.

Available in 4 sizes : 76 x 45 cm, 76 x 90 cm, 76 x 120 cm and 76 x 150 cm.  
Max load 300 kg .



L4a. SC ramp in use



L4b. SC ramp folded for transport

**L5. Double-folded ramps - DF-ramps**

Made from aluminium and covered with a friction surface. Consists of 2 ramps that are hinged at the middle and carried like a suitcase.

Available in 3 sizes, 76 x 180 cm, 76 x 240 cm and 76 x 300 cm. Max. load 300 kg.



**L5a.** DF-ramp in use



**L5b.** DF-ramps folded for transport

**L6. Suitcase ramps - SR ramps**

Made from aluminium with a friction surface. Fourfold foldable one-piece ramp, complete with carrier handle.

Available in 2 sizes: 76 x 180 cm and 76 x 240 cm .  
Max. load 300 kg.



**L6a.** SR ramp in use



**L6b.** The SR ramp is fourfold foldable for transport

**L7. Doorstep ramp - TSR Ramp**

Made of aluminium with anti-slip pattern. Suitable for terraces, entrances or in the garden. It can be height adjusted very precisely from 15 to 22.5 cm.

There are 2 types: 80 x 100 cm with 30 cm flap, and 80 x 100 cm with 40 cm flap.  
Max. load 225 kg.



**L7a.** TSR Ramp in use



**L7b.** TSR Ramp: Height adjustment

## L8. Terrace ramps



**L8a.** Terrace ramp



**L8b.** Height adjustment mechanism

Terrace ramp made of aluminium with anti-slip pattern. The height is easily adjusted by loosening a bolt on each side. The ramp is designed so that it is placed on the ground - a permanent installation is not necessary. It is mounted so that it reaches the upper level of the threshold. Since it has no side-protection it is also possible to open the door while the ramp is in use.

The ramp is available in 3 sizes (see below).  
Max. load 300 kg.

	<b>Size (W x L)</b>	<b>Height adjustment</b>
<b>Terrace ramp T80060</b>	80 x 60 cm	5 - 9 cm
<b>Terrace ramp T80075</b>	80 x 75 cm	8 - 13 cm
<b>Terrace ramp T80125</b>	80 x 125 cm	11 - 28 cm