



Foldable shower-chair

Frame made from white coloured aluminium with ergonomic seat with cut-away for better personal hygiene.

OX1. Foldable shower chair

Measurements:...52,5 x 50 x 80-90 cm (W x d x H)
Folded size:.....55,5 x 16,5 x 80,5 cm (W x D x H)
Seat size:.....44 x 39 cm (W x D) with cut-away of 12,5 x 18 cm (W x D)
Height:.....adjustable 44,5 - 54,5 cm in 5 positions
Armrest:.....height 22 cm
Backrest:.....12 x 36 cm (W x D), adjustable up to a height of 33.5 cm above the seat
Weight:.....6,4 kg
Max. load:.....100 kg



Foldable shower-chair with WC

Frame made from white coloured aluminium with a removable toilet bucket with lid and carrying handle.

OX2. Foldable shower-chair with WC

Measurements:...52,5 x 50 x 80-90 cm (W x D x H)
Folded size:.....55,5 x 16,5 x 80,5 cm (W x D x H)
Seat size:.....30 x 50 (W x D)
Height:.....adjustable from 44,5 to 54,5 cm in 5 positions
Armrests:.....height 22 cm
Backrest:.....12 x 36 cm (W x D), adjustable up to a height of 33.5 cm above the seat
Weight:.....7,4 kg
Max. load:.....100 kg



OX3. Shower chair, round Ø 33 cm

Aluminium shower chair with round white plastic seat adjustable height.

Size:.....seat Ø 33 cm
Width:.....33 cm
Height:.....adjustable 40-50 cm
Weight:.....1,66 kg
Max. load:.....130 kg



OX4. Wall mounted shower-chair

Wall mounted shower chair, cut-away for better personal hygiene.

Size:.....44 x 39 cm (W x D)
Cut-away of:.....12,5 x 18 cm (W x D)
Height:.....adjustable 46-53,5 cm
Weight:.....3,5 kg
Max. load:.....100 kg



OX5. Bath board

Made from injection moulded white plastic. Anatomic formed with non-slip surface and handle.

Size (L x W x H):.....70,5 x 32 x 18 cm
For bathtubs:.....41-64 cm (inside)
Width adjustment:...40,5 - 65 cm
Weight:.....2,7 kg
Max. load:.....190 kg



OX6. Sliding bathtub seat

Sliding bathtub seat

Bathtub seat to use for easy entry and exit from the bathtub. Aluminium frame with seat, back and soap dish made of white plastic. The seat is adjustable in height.

Total width:.....106 cm
Depth:.....46,5 cm
Height:.....adjustable 79,5 - 89,5 cm
Seat-width:.....50 cm
Seat-depth:.....39,5 cm
Seat-height:.....adjustable 46 - 56 cm
Weight:.....6,50 kg
Max. user weight:..130 kg



OX7. Shower chair, grey, with soft seat

Shower chair made of grey metallic coloured aluminium with a soft seat of grey PU. The seat is adjustable in height. The chair is also available with the frame in white coloured aluminium.



OX8. Back for bath-chair OX7 - accessory

Available as accessory is a soft backrest mountable without tools. Available in frame colours white or grey.

	Seat width 42 cm	Seat width 47 cm	Seat width 52 cm
Size (WxDxH)	52,5 x 53,5 x 60-77 cm	57,5 x 53,5 x 60-77 cm	62,5 x 53,5 x 60-77 cm
Seat size (WxD)	42 x 34 cm	47 x 34 cm	52 x 34 cm
Seat height adjustable	42 - 59 cm	42 - 59 cm	42 - 59 cm
Weight	2,70 kg	2,90 kg	3,30 kg
Max. Load	150 kg	150 kg	150 kg



OX9. Bathtub seat - grey with soft polyurethane seat

Bathtub seat with hygiene cut-away. Silver-coloured aluminium frame. Seat and backrest (optional) of soft grey polyurethane. The optional backrest is fixed with clips and is easily removable.

Size:.....76-91 x 42 x 23 cm (W x D x H)
Seat size:.....41 x 34 cm with cut-away
Height:.....backrest: 38 cm
For bathtubs:.....from 57 to 66 cm (inner width)
Weight:.....2,6 kg
Max. Load:.....150 kg



OX10. Swivel bather - grey with soft polyurethane seat

Swivel bather with hygiene cut-away. Silver-coloured aluminium frame. Seat and backrest are made of soft grey PUR. The backrest and the armrests are fixed with clips and are easily removable.

Size:.....71 x 54 x 57 cm (W x D x H)
Seat size:.....41 x 34 cm with cut-away
Height:.....armrests: 20 cm - backrest: 38 cm
For bathtubs:.....from 57 to 66 cm (inner width)
Weight:.....4,88 kg
Max. Load:.....130 kg