

Standard seat cushions

Standard wheelchair cushion made from high quality cold foam with an especially high load capacity, or with visco-elastic foam. The lower part of the cushion is covered with an anti-slip material. In this way the cushion can be used in any wheelchair without the need to fit with existing velcro strips. The cover is made from Polyester, durable, washable and elastic, and is fitted with a zipper.

The cushion is 5 cm high and all have a seat depth of 43 cm.
Deliverable in 6 widths from 37,5 cm to 50 cm.

Seat cushion with PU foam or high density visco-elastic foam and black polyester cover



P1. The cushion seen from the front.



P2. Anti-slip material



Seat cushion with PU foam or high density visco-elastic foam and black incontinence cover



P3. Cushion with incontinence cover

The cushion is also available with incontinence cover

Note:

Available also with individual imprint for an extra charge.

Seat cushion with a double layer (PU foam and visco-elastic foam)



P4a. Incontinence cover



P4b. Anti-slip material



P4c. Viscoelastic foam

Seat cushion made from one layer of normal PU foam and one layer of visco-elastic foam. It comes with an incontinence cover.

The cushion is deliverable in 5 different sizes (W x D x H):
40.0 X 40 X 8 cm, 42.5 X 43 X 8 cm, 45.0 X 45 X 8 cm, 47.5 X 45 X 9 cm and 50.0 X 48 X 9 cm

**P5.** Wheelchair cushion with valve for regulation of air pressure**Mobilex air cushion**

The Mobilex air-cushion is equipped with a special valve for air pressure. The user "immerses" in the foam, which conforms to the user and distributes the user's weight over the surface of the cushion. By closing the valve the foam maintains its form for an optimal pressure distribution.

Sizes (W x D x H):

40 x 40 x 8 cm

45 x 45 x 8 cm

50 x 48 x 8 cm

**P6.** Swivel seat with Ø 45 cm**Mobilex swivel seat**

Can be used in a car or other place where rotatability is wanted.

Diameter : 45 cm

Max. load : 100 kg

**P7.** Wedge with incontinence cover**Mobilex wedge with Decutex cover**

Wedge in visco foam with Decutex cover

Size: (B x T x H): 45 x 43 x 10/0 cm

**P8.** Wedge made of regenerated foam**Standard wedge**

Made of regenerated foam

Size (B x D x H):

45 X 45 X 8/1 cm

45 X 45 X 10/1 cm