

The "Gazelle" lightweight rollator has a frame made of silver coloured aluminium and can be folded along the direction of travel, such as a wheelchair. Standard equipped with a basket and a comfortable backrest. The large wheels with a diameter of 200 mm ensure stability and good rolling ability. The seat and the handles can be adjusted in height. He has efficient brakes that also can be blocked and is secured against unwanted folding.

Technical data	Gazelle Standard item nr. 312007	Gazelle small item nr. 312006	Gazelle large item nr. 312008
Hand grip height min./max.	83-104 cm	76-97 cm	83-104 cm
Seat height min./max.	55-60 cm	48-53 cm	55-60 cm
Seat width	46 cm	46 cm	51 cm
Outer width	56 cm	56 cm	61 cm
Weight	7,6 Kg	6,9 Kg	7,8 Kg
Max. load	120 Kg	120 Kg	120 Kg



RX1a. Rollator "Gazelle" Standard



RX1b. Seat-height adjustment



RX1c. The "Gazelle" rollator folded



R3. Speed reducing device



R4. Brake handle



R5. Brake handle for one-arm drive with cables



R6. Clip to fix cables to the rollator



R7. Brake shoe



R8. Crutch holder