

Further in the Mobilex program

A large selection of mobile aluminium wheelchair ramps.



Telescopic ramps, single or double folded ramps, treshold ramps and terrace ramps.

Available in different sizes.

Ask your dealer



The company Mobilex A/S was founded in 2003 and we offer a wide range of aids for handicapped people like wheelchairs, walking aids, ramps, bath products and a wide range of spare parts.

We distribute our products throughout Europe and serve our customers through our own sales representatives.


High quality standards, modern technology, fast delivery, a customer-oriented service, and motivated employees are the key of our success.

We look forward to a pleasant cooperation.



Mobilex A/S
Grønlandsvej 5
8660 Skanderborg
Denmark

Tel. +45 8793 2220
Fax +45 8793 1777
info@mobilex.dk
www.mobilex.dk



RUBBER RAMPS




Available in different heights from 4 to 60 mm and in widths of 700, 900 and 1000 mm



Available sizes

The Mobilex rubber ramps are made of rubber (TPE) and contain no PVC. They are available in different heights (see table). The standard width is 900 mm. Also available in widths of 700 or 1000 mm.



The Mobilex rubber ramps are fitted with double-sided adhesive tape on the underside for attachment to smooth floors. For laying on carpets or similar rubber ramps can also be supplied without adhesive tape.

Rubber ramps 700 mm

Item nr.	Description (H x L x W)	Adhesive
30610L	Rubber ramp 10 x 100 x 700 mm	yes
30612L	Rubber ramp 12 x 120 x 700 mm	yes
30616L	Rubber ramp 16 x 150 x 700 mm	yes
30620L	Rubber ramp 20 x 150 x 700 mm	yes
30624L	Rubber ramp 24 x 150 x 700 mm	yes
30628L	Rubber ramp 28 x 150 x 700 mm	yes
30632L	Rubber ramp 32 x 200 x 700 mm	yes
30636L	Rubber ramp 36 x 200 x 700 mm	yes
30640L	Rubber ramp 40 x 250 x 700 mm	yes
30644L	Rubber ramp 44 x 250 x 700 mm	yes
30648L	Rubber ramp 48 x 250 x 700 mm	yes

Rubber ramps 900 mm

Item nr.	Description (H x L x W)	Adhesive
306040	Rubber ramp 4 x 40 x 900 mm	yes
306060	Rubber ramp 6 x 60 x 900 mm	yes
306080	Rubber ramp 8 x 80 x 900 mm	yes
306100	Rubber ramp 10 x 100 x 900 mm	yes
306120	Rubber ramp 12 x 120 x 900 mm	yes
306160	Rubber ramp 16 x 150 x 900 mm	yes
306200	Rubber ramp 20 x 150 x 900 mm	yes
306240	Rubber ramp 24 x 150 x 900 mm	yes
306280	Rubber ramp 28 x 200 x 900 mm	yes
306320	Rubber ramp 32 x 200 x 900 mm	yes
306360	Rubber ramp 36 x 200 x 900 mm	yes
306400	Rubber ramp 40 x 250 x 900 mm	yes
306440	Rubber ramp 44 x 250 x 900 mm	yes
306480	Rubber ramp 48 x 250 x 900 mm	yes
306520	Rubber ramp 52 x 300 x 900 mm	yes
306560	Rubber ramp 56 x 300 x 900 mm	yes
306600	Rubber ramp 60 x 300 x 900 mm	yes

306041	Rubber ramp 4 x 40 x 900 mm	no
306061	Rubber ramp 6 x 60 x 900 mm	no
306081	Rubber ramp 8 x 80 x 900 mm	no
306101	Rubber ramp 10 x 100 x 900 mm	no
306121	Rubber ramp 12 x 120 x 900 mm	no
306161	Rubber ramp 16 x 150 x 900 mm	no
306201	Rubber ramp 20 x 150 x 900 mm	no
306241	Rubber ramp 24 x 150 x 900 mm	no
306281	Rubber ramp 28 x 200 x 900 mm	no
306321	Rubber ramp 32 x 200 x 900 mm	no
306361	Rubber ramp 36 x 200 x 900 mm	no
306401	Rubber ramp 40 x 250 x 900 mm	no
306441	Rubber ramp 44 x 250 x 900 mm	no
306481	Rubber ramp 48 x 250 x 900 mm	no
306521	Rubber ramp 52 x 300 x 900 mm	no
306561	Rubber ramp 56 x 300 x 900 mm	no
306601	Rubber ramp 60 x 300 x 900 mm	no


Rubber ramps 1000 mm

Item nr.	Description (H x L x W)	Adhesive
306042	Rubber ramp 4 x 40 x 1000 mm	yes
306062	Rubber ramp 6 x 60 x 1000 mm	yes
306082	Rubber ramp 8 x 80 x 1000 mm	yes
306102	Rubber ramp 10 x 100 x 1000 mm	yes
306122	Rubber ramp 12 x 120 x 1000 mm	yes
306162	Rubber ramp 16 x 150 x 1000 mm	yes
306202	Rubber ramp 20 x 150 x 1000 mm	yes
306242	Rubber ramp 24 x 150 x 1000 mm	yes
306283	Rubber ramp 28 x 200 x 1000 mm	yes

306282	Rubber ramp 28 x 200 x 1000 mm	no
306322	Rubber ramp 32 x 200 x 1000 mm	no
306362	Rubber ramp 36 x 200 x 1000 mm	no
306402	Rubber ramp 40 x 250 x 1000 mm	no
306442	Rubber ramp 44 x 250 x 1000 mm	no
306482	Rubber ramp 48 x 250 x 1000 mm	no
306522	Rubber ramp 52 x 300 x 1000 mm	no
306562	Rubber ramp 56 x 300 x 1000 mm	no
306602	Rubber ramp 60 x 300 x 1000 mm	no

Accessories

For mounting



The fixing set for rubber ramps consists of two holding points with which the rubber ramps (from 28 mm height) can be attached to the door frame.

For adaptation

Saw blade for shortening



Saw blade for standard jigsaws to adjust the width of Mobilex rubber ramps. The Mobilex rubber ramps can easily be shortened.

Double adhesive tape



Mobilex rubber ramps are supplied with adhesive strips on the underside as standard. A double adhesive tape is available for replacement or conversion.

Solvent for adhesive removal



A special solvent is available as an accessory to remove the double-sided adhesive tapes from the underside of the rubber ramps.