

THRESHOLD RAMPS

DS Ramps

Aluminium doorstep ramps with ribbed surface. The DS2176 has three adjustable feet. The DS4076 is foldable and has 2 milled handles and a hinge for folding.




DS2176:
Length: 76 cm
Width: 21 cm
Weight: 1,55 kg
Max. load. 350 kg
Max. elevation/height: 3-7 cm

DS4076:
Length: 76 cm
Width: 40 cm
Weight: 3,3 kg
Max. load 350 kg
Max elevation/height 3-5 cm

DS double ramp


The DS4025 threshold ramp is a set of 2 hinged parts made from aluminium with an anti-slip surface. The ramp comes with a useful bag for transport.



DS4025:
2 pcs. 2,5 x 40 x 25 cm (H x D x W)

Rubber Ramps

Rubber ramps with anti-slip surface for use among others in bathrooms and at doorsteps. The ramps come in 2 versions: with glue protected by a folio and without glue. The ramps are available in widths of 700mm, 900 mm or 1000 mm. If the ramps in standard measure are too long they can be easily cut with a Stanley knife.




Available in various heights from 4 to 60 mm and widths from 700, 900 or 1000 mm
Please ask for the info sheet for the rubber ramps

TELESCOPIC RAMPS

TM Ramps

The ramps are made of aluminum. The perforated lanes provide an ideal balance between weight and durability. In addition, the permeability is advantageous in rain, dirt or snow. The height of the side profile provides sufficient grip and directional stability. The ramps consist of two separate drive-on rails. The model TM-055 consists of one single segment. The model TM-150 and TM-210 consists of 2 telescopic segments, model TM-300 has also a further element which is foldable. In the folded state, the ramps are secured by a clip closure. The TM ramps must be used in pairs.



TM Ramps	TM-055	TM-150	TM-210	TM-300
Length	55 cm	150 cm	210 cm	300 cm
Width (per rail)	25 cm	25 cm	25 cm	26 cm
Surface width	20 cm	20 cm	20 cm	20 cm
Height of side protect.	5,5 cm	5,5 cm	5,5 cm	5,5 cm
Max. Elevation / height	12° / 11 cm	12° / 30 cm	12° / 40 cm	12° / 60 cm
Weight	2 x 1,45 kg	2 x 3,95 kg	2 x 5,8 kg	2 x 8,65 kg
Max. load (one pair)	300 kg	300 kg	300 kg	250 kg
Folded size cm (LxWxH)	55 x 25 x 5,5	98 x 27x 9	128 x 27 x 9	108 x 26 x 13,5

TR/LR Ramps


The ramps are made of aluminium with a friction surface. They consist of two separate drive-on rails with 3 segments which are provided with a locking mechanism to prevent unintentional extension and retraction of the segments. The models "LR" are equipped with a convenient carrying case that can be attached to the wheelchair. The carrying case is also available for the TR model as an optional accessory (surcharge). The telescopic ramps must be used in pairs.



TR/LR Ramps	TR-150	TR-210	TR-300	LR-150	LR-210
Length	150 cm	210 cm	300 cm	150 cm	210 cm
Width (per rail)	26 cm	26 cm	26 cm	16 cm	16 cm
Surface width	19 cm	19 cm	19 cm	11 cm	11 cm
Height of side protect.	7/5 cm	7/5 cm	7/5 cm	5,5/3 cm	5,5/3 cm
Max. Elev. / height	12° / 30 cm	12° / 40 cm	12° / 60 cm	12° / 30 cm	12° / 40 cm
Weight	2 x 5,5 kg	2 x 7,5 kg	2 x 9,5 kg	2 x 3 kg	2 x 4 kg
Max. load (one pair)	225 kg	225 kg	225 kg	180 kg	180 kg
Folded size (LxWxH)	74 x 26,5 x 7	95x26,5x7	125x26,5x7	74x16,5x5,5	95x6,5x5,5



WHEELCHAIR RAMPS



Dealer

SINGLE FOLDED RAMPS

SINGLE FOLDED RAMPS

SC Ramps

Inexpensive suitcase ramp in aluminium with anti slip surface. Hinged at the middle and carried like a suitcase.



SC Ramps	SC-045	SC-090	SC-120	SC-150	SC-180
Length	46 cm	92 cm	122 cm	152 cm	182 cm
Width	77 cm	77 cm	77 cm	77 cm	77 cm
Surface width	76 cm	76 cm	76 cm	76 cm	76 cm
Height of side protect	3,5 cm	3,5 cm	3,5 cm	3,5 cm	3,5 cm
Max. elevation / height	12° / 10 cm	12° / 20 cm	12° / 25 cm	12° / 30 cm	12° / 36 cm
Weight	4,50 kg	9,00 kg	12,00 kg	15,00 kg	18,00 kg
Max. load	300 kg	300 kg	300 kg	300 kg	300 kg
Folded size (L x W x H)	46 x 40 x 8	92 x 40 x 8	122 x 40 x 8	152 x 40 x 8	182 x 40 x 8

SF Ramps

The ramps are made of aluminium with a friction surface and are collapsible in the longitudinal direction - a so-called flat-ramp. This offers the opportunity for volunteers to be able to stand on the ramp together with the wheelchair and the patient. The ramps are equipped with a carrying handle and hinged at the middle so they are portable like a suitcase.




SF Ramps	SF-040	SF-090	SF-120	SF-150
Length	45 cm	90 cm	120 cm	150 cm
Width	71 cm	80 cm	80 cm	80 cm
Surface width	71 cm	76 cm	76 cm	76 cm
Height of side protect.	--	6 cm	6 cm	6 cm
Max. Elevation / height	12° / 10 cm	12° / 20 cm	12° / 25 cm	12° / 30 cm
Weigh	4 kg	8 kg	10,5 kg	12,5 kg
Max. load	300 kg	300 kg	300 kg	300 kg
Folded size (L x W x H)	45 x 37 x 5	90 x 40 x 10	120 x 40 x 10	150 x 40 x 10

DOUBLE FOLDED RAMPS

DOUBLE FOLDED RAMPS

DF Ramps


Made from aluminium and covered with a friction surface. Consists of 2 ramps that are hinged at the middle and can be carried like a suitcase.



DF-Ramps	DF-180	DF-240	DF-300
Length	180 cm	240 cm	300 cm
Width	79 cm	79 cm	79 cm
Surface width	75 cm	75 cm	75 cm
Height of side protect.	6 cm (inside)	6 cm (inside)	6 cm (inside)
Max. Elevation / height	12° / 36 cm	12° / 50 cm	12° / 60 cm
Weight	16 kg	22 kg	27,2 kg
Max. load	300 kg	300 kg	300 kg
Folded size (L x W x H)	91 x 36 x 19	124 x 40 x 19	154 x 40 x 19

SR Ramps

Made from aluminium with a friction surface. Fourfold foldable one-piece ramp, complete with carrier handle.



SR-Ramps	SR-180	SR-240
Length	180 cm	240 cm
Width	80 cm	80 cm
Surface width	72 cm	72 cm
Height of side protect.	Outer side: 6 cm	Outer side: 6 cm
Max. Elevation / height	12° / 36 cm	12° / 50 cm
Weight	14,1 kg	18,5 kg
Max. load	300 kg	300 kg
Folded size (L x W x H)	102 x 28 x 38	135 x 31 x 38

CHANNEL RAMP

CHANNEL RAMP

CH Ramp

The ramp is made of aluminium with a friction surface and consists of two separate drive-on rails that are hinged at the middle. Folded the ramp can be carried with 2 handles like a suitcase. The CH ramps must be used in pairs.




CH-180:
Length: 180 cm
Rail width: 23,5 cm
Surface width: 21 cm
Weight: 13 kg (1 pair)
Max. load: 300 kg (1 pair)
Max. elevation/height: 35 cm

TERRACE RAMPS

TERRACE RAMPS

TSR Ramps


Made of aluminium with anti-slip pattern. Suitable for terraces, entrances or in the garden. It can be height adjusted very precisely.



TSR-130/TRS-140:
Length: 130 or 140 cm
Width: 80 cm
Width of the flap: 75 cm
Weight: 14,7 or 15,6 kg
Max. load 225 kg
Max. elevation/height 20 cm

T Ramps

Terrace ramp made of aluminium with anti-slip pattern. The height is easily adjusted by loosening a bolt on each side. It is mounted so that it reaches the upper level of the threshold. Since it has no side-protection it is also possible to open the door while the ramp is in use.



T80060/T80075/T80080:
Max. load 225 kg
Height: 5-9, 8-13 or 11-28 cm
Weight: 8 kg, 12 kg or 21,5 kg
Max. elevation/height: 9, 13 or 28 cm
Length x Width: 60, 75 or 125 cm x 80 cm