

The new "Bison" verticalizer and walking frame in modern design



Technical data

Width: 60-80 cm
Height: 86-128 cm
Length: 81 cm
Weight: 16,4 kg
Max. load: 150 kg

The "Bison" walking frame is a tool for disabled or elderly people who need additional support to maintain balance or stability while walking. The frame is made of powder coated aluminium and offers various adjustment functions: the armrests are supported by a gas spring to grant easy height adjustment. The arm rests are adjustable laterally. The angle and depth of the handles are also adjustable. All settings can be done without tools.

