

Indoor rollator



312410



312411



312412



312413



312415 - SUPREME

- Aluminium frame which folds easily for storage.
- Lightweight rollator that is very stable and easy to manoeuvre.
- Equipped with four smooth-running wheels.
- Equipped with two lockable brakes, which are operated by means of a handle.
- The rollator is tested and approved for a user weight of up to 100 kg.
- The rollator is only approved for indoor use.

As standard the Jaguar comes with:

- a bag
- a removable tray

Rollator Jaguar	Various colours	SUPREME
Height of push handles	83,5 - 94,5 cm	90,3 - 96 cm
Folded rollator size (H x W x D)	83,5 x 54,5 x 27 cm	90,3 x 56,5 x 27 cm
Width	54,5 cm	56,5 cm
Length opened / folded	66 cm / 27 cm	66 cm / 27 cm
Weight	6,3 kg	6,0 kg
Max. load (rollator / bag)	100 kg / 5 kg	100 kg / 5 kg
Brakes	2 lockable with parking function	
Wheel sizes	Front and back: 7" (175 x 27 mm)	



312414



312416



312417