

TIGER

OUTDOOR ROLLATOR

The **Mobilex Arthritis Rollator Tiger Outdoor** is based on the elegant Rollator Lion and therefore has all its features and functions. As an Arthritis Rollator it also has soft and comfortable forearm supports. These can be adjusted in depth, angle and height and can be turned inwards in 2 steps of 22.5 degrees each. The rollator is extremely stable and easy to manoeuvre. It has large, smooth-running soft wheels and easy to operate and lockable brakes. In driving and transport position it is secured against unintentional unlocking. A detachable shopping bag and a special back strap are included.

The Arthritis Rollator Tiger Outdoor is also approved for **outdoor use** and for a user weight of up to **200 kg**.




TIGER

OUTDOOR ROLLATOR

The Mobilex Arthritis Rollator Tiger Outdoor is available in two colours, **metallic brown** and **anthracite grey**.

The rollator has also been **tested and approved for outdoor use** and for a user weight of **200 kg**.



On the Arthritis Rollator Tiger Outdoor, the **forearm supports can be rotated** inwards in 2 steps of 22.5 degrees each.

Mobilex Rollator Tiger Outdoor	Antracite grey Item nr. 312205	Brown metallic Item nr. 312206
<i>Height of the forearm support (adjustable)</i>	87 - 110 cm	87 - 110 cm
<i>Seat height</i>	62 cm	62 cm
<i>Total width opened / folded</i>	71 cm / 31 cm	71 cm / 31 cm
<i>Inner width between push handles</i>	56 cm	56 cm
<i>Length (depth)</i>	77 cm	77 cm
<i>Weight (+ bag 380 g)</i>	12 kg	12 kg
<i>Max. load (rollator / bag)</i>	200 kg / 5 kg	200 kg / 5 kg
<i>Brakes</i>	2 lockable with parking function	
<i>Wheel sizes</i>	Front and back 200 x 35 mm	