

# ROLLATORS

MOBILEX SELECTION OF ROLLATORS





# TABLE OF CONTENTS

LEOPARD rollator.....	3
GEPARD rollator.....	4
BUFFALO HD rollator.....	4
ROLLOGUARD.....	5
PANDA COMFORT rollator.....	5
TIGER & TIGER OUTDOOR rollator.....	6
LION rollator.....	7
PANTHER rollator.....	7
PUMA AIR rollator.....	8
JAGUAR rollator.....	9
IMPALA rollator.....	9
ACCESSORIES.....	10
BISON walking frame.....	11
DYNAWALK walking frame.....	11

---

## General description

Mobilex rollators are of high quality, user-friendly, stable, easy-to-use and contribute to “a better life in motion” for the users. We have the right rollator for every situation and challenge. Whatever you’re looking for, you’ll find it here. Our extensive range of rollators meet the diverse needs of users.

Although our rollators are named after speedy animals, its up to the user to decide the speed as well as the model and colour. Naturally, all Mobilex rollators are equipped with well-functioning and user-friendly brakes for the users safety and comfort. A lightweight and very stable frame construction provides the necessary stability during use. Mobilex rollators are available up to a user weight of 200 kg and a seat width of up to 56 cm. With a maximum height of the drive handles up to 112 cm, even the tallest user will be able to find a suitable model. Whether for indoor or outdoor use, we at Mobilex A/S offer rollators for all needs. We also have several modern and elegant models for patients with arthritis in our wide repertoire.

On the following pages you can immerse yourself in our wide range of Mobilex rollators and the extensive accessory program. If you do not find what you are looking for, you are very welcome to contact us or one of our many qualified Mobilex dealers.

### More information?

On these pages you will find a small selection of our extensive product range. If you are interested in our full range, we invite you to visit our website [www.mobilex.dk](http://www.mobilex.dk), where you can view all products and download our catalog. The QR code will take you to our website.



# LEOPARD ROLLATOR



The Mobilex lightweight rollator **Leopard** has an aluminium frame that can be easily folded sideways, which in folded state reduces the required space to a minimum. The rollator is equipped with a crutch holder, a back belt and a handy removable bag as standard. The bag has carrying handles and can be used as a shopping bag. Easy to open and fold, the rollator is suitable for use both outdoors and indoors. The comfortable seat is made of tear resistant nylon. The push handles can easily be adjusted in height without use of tools. The rollator comes with two lockable brakes on the rear wheels and is tested and approved for a user weight of up to 150 kg.



Rollator Leopard	Standard 312301 (red) 312304 (anthracite) 312309 (blue) 312300 (green) 312303 (platin)	Standard+Softwheels 312308 (red) 312306 (anthracite) 312314 (blue) 312313 (green) 31230A (platin)	Small 312302 (anthracite) 312319 (red)	X-Small 31230K (anthrazit)	Small and Softwheels 312318 (anthracite) 31230E (red)
Height of the push handles (adjustable)	74 - 102 cm	74 - 102 cm	66 - 86 cm	63 - 79 cm	66 - 86 cm
Recommended user heights	150 - 200 cm	150 - 200 cm	135 - 170 cm	130 - 165 cm	135 - 170 cm
Seat height	62 cm	62 cm	55 cm	50 cm	55 cm
Total width opened / folded	60 / 22,5 cm	60 / 22,5 cm	60 / 22,5 cm	60 / 22,5 cm	60 / 22,5 cm
Inner width between push handles	46 cm	46 cm	46 cm	46 cm	46 cm
Length (depth)	66 cm	66 cm	66 cm	66 cm	66 cm
Weight (+ bag 380 g)	7,4 kg	7,6 kg	7,1 kg	6,9 kg	7,3 kg
Max. load (rollator / bag)	150 kg / 5 kg	150 kg / 5 kg	150 kg / 5 kg	150 kg / 5 kg	150 kg / 5 kg
Brakes	2 lockable with parking function				
Wheel sizes	Front and back 200 x 35 mm				

# GEPARD ROLLATOR

The Mobilex lightweight rollator **Gepard** has a frame made of carbon fibre. It can be folded sideways, which in folded state reduces the required space to a minimum. The rollator comes with a crutch holder, a comfortable soft seat and a removable bag. The rollator is easy to open and fold and with its large wheels, it is suitable for use both outdoors as indoors. The push handles can easily be adjusted in height. The rollator comes with two lockable brakes on the rear wheels and is tested and approved for a user weight of up to 150 kg (200 kg for the HD model).



Gepard Rollator	Standard 312401 (black) 312403 (yellow) 312404 (purple/black)	With SOFT wheels 312406 (black) 312407 (yellow) 312408 (purple/black)	Small 312402 (black)	HD * 312405 (black)
Height of the push handles (adjustable)	74 - 102 cm	74 - 102 cm	66 - 86 cm	74 - 102 cm
Recommended user heights	150 - 200 cm	150 - 200 cm	135 - 170 cm	150 - 200 cm
Seat height	62 cm	62 cm	55 cm	62 cm
Total width opened / folded	60 / 22,5 cm	60 / 22,5 cm	60 / 22,5 cm	68 / 22,5 cm
Inner width between push handles	46 cm	46 cm	46 cm	56 cm
Length (depth)	66 cm	66 cm	66 cm	66 cm
Weight (+ bag 380 g)	5,4 kg	5,5 kg	5,1 kg	5,7 kg
Max. load (rollator / bag)	150 kg / 5 kg	150 kg / 5 kg	150 kg / 5 kg	200 kg / 5 kg
Brakes	2 lockable with parking function			
Wheel sizes	Front and back: 200 x 35 mm			

# BUFFALO HD ROLLATOR

The Mobilex lightweight rollator **Buffalo HD** has an aluminium frame that can be folded sideways, which in folded state reduces the required space to a minimum. Its structure is reinforced with additional elements. The rollator is equipped with a back belt, a crutch holder and a handy removable bag. The bag has carrying handles and can also be used as a shopping bag. The rollator is easy to open and fold and suitable for use both outdoors and indoors. The comfortable seat is made of tear resistant nylon. The push handles can easily be adjusted in height without use of any tools. The rollator comes with two lockable brakes on the rear wheels and is tested and approved for a user weight of up to 200 kg.



Buffalo HD Rollator	Standard 312305	Small 312307	Standard with SOFT wheels 312385	Small with SOFT wheels 312387
Height of the push handles (adjustable)	74 - 102 cm	66 - 86 cm	74 - 102 cm	66 - 86 cm
Recommended user heights	150 - 200 cm	125 - 170 cm	150 - 200 cm	125 - 170 cm
Seat height	62 cm	55 cm	62 cm	55 cm
Total width opened / folded	68 / 22,5 cm	68 / 22,5 cm	68 / 22,5 cm	68 / 22,5 cm
Inner width between push handles	56 cm	56 cm	56 cm	56 cm
Length (depth)	66 cm	66 cm	66 cm	66 cm
Weight (+ bag 380 g)	7,8 kg	7,5 kg	7,9 kg	7,6 kg
Max. load (rollator / bag)	200 kg / 5 kg	200 kg / 5 kg	200 kg / 5 kg	200 kg / 5 kg
Brakes	2 lockable with parking function			
Wheel sizes	Front and back: 200 x 35 mm			



# ROLLOGUARD

The innovative and fall-preventive braking systems from RolloGuard are perfectly combined with our high-quality Mobilex rollators. RolloGuard's fall prevention for rollators reduces the risk of crashes or fall injuries with a rollator. The patented fall prevention technology prevents accidental movements and thus enables safer mobility. RolloGuard's Instant Fall Prevention is **compatible with the Buffalo, Gepard, Leopard and Tiger rollator**. The RolloGuard fall prevention uses the "Transversor" technology, which protects against unintentional movements and makes mobility safer - for a better and safer life in motion.

## Adjustable protection level

A rollator equipped with RolloGuard is marked with a sticker on the brake drum. This sticker is referenced with a specific color, description and intended use. A distinction is made between the three support levels; high (red), medium (yellow) and low (green), all of which have different sensitivities, which affects the activation of the fall protection.



# PANDA COMFORT ROLLATOR

The lightweight rollator **Panda Comfort** has a dark blue aluminium frame. The rollator folds easily sideways and when folded it takes up minimal space, making it easy to transport. As standard, the rollator is equipped with a back belt, which provides safety and support when the seat is used. The seat is made of heavy-duty, durable nylon. A closed shopping bag is fitted to the front, which has two handles so that it can easily be removed and carried separately. The large soft wheels make the rollator good and stable for both indoor and outdoor use. In addition, the large wheels also provide good shock absorption. The push handles can easily be adjusted in height without the use of any tools. The rollator is equipped with two lockable brakes that brake on the rear wheels. Two integrated tilting pedals are fitted near to the rear wheels to facilitate overcoming of minor obstacles. As standard the Panda Comfort comes with a transport lock, a removable shopping bag, a back belt and a crutch holder. Otherwise, a wide range of other accessories is available. The rollator is tested and approved for a user weight of up to 150 kg.



<b>Panda Comfort Rollator</b>	<b>312601</b>
<i>Height of the push handles (adjustable)</i>	77,5 - 99,1 cm
<i>Recommended user heights</i>	160 - 200 cm
<i>Seat height</i>	59 cm
<i>Total width opened / folded</i>	65 cm / 26,4 cm
<i>Inner width between push handles</i>	45,5 cm
<i>Length (depth)</i>	77,5 cm
<i>Weight</i>	9,44 kg
<i>Max. load (rollator / bag)</i>	150 kg / 5 kg
<i>Brakes</i>	2 lockable with parking function
<i>Wheel sizes</i>	Front and back: 285 mm x 45 mm



# TIGER / TIGER OUTDOOR ROLLATOR



The Mobilex Arthritis **Tiger** and **Tiger Outdoor** rollators are based on the elegant Leopard (Tiger) and Lion (Tiger Outdoor) rollators, respectively, and therefore also have their characteristics and functions. In addition, the soft and comfortable forearm supports with their brake handles are adjustable in depth, angle and height. Apart from the model 312201, the armrests can also be turned inwards in two angles by 22.5 degrees each. The rollators are extremely stable and easy to manoeuvre. They have large, smooth-running wheels and lockable brakes. In driving position they are secured against unintentional unlocking. As standard, the rollators have a transport lock and a removable shopping bag. A wide range of accessories is also available. The rollator Tiger Standard is only approved for indoor use on level floors and for a user weight of up to 150 kg. The Tiger Outdoor model is also approved for outdoor use with a user weight of up to 200 kg.

<b>Tiger Arthritis Rollator</b>	<b>Standard</b> 312201 (anthracite grey) 312211 (red)	<b>With rotatable armrests</b> 312202 (metallic brown)
<i>Height of the push handles (adjustable)</i>	90 cm - 115 cm	90 cm - 115 cm
<i>Recommended user heights</i>	150 - 200 cm	150 - 200 cm
<i>Seat height</i>	62 cm	62 cm
<i>Total width opened / folded</i>	62 / 22,5 cm	62 / 22,5 cm
<i>Inner width between push handles</i>	46 cm	46 cm
<i>Length (depth)</i>	66 cm	66 cm
<i>Weight (+ bag 380 g)</i>	9,5 kg	9,5 kg
<i>Max. load (rollator / bag)</i>	150 kg / 5 kg	150 kg / 5 kg
<i>Brakes</i>	2 lockable with parking function	
<i>Wheel sizes</i>	Front and back: 200 x 35 mm	



<b>Tiger Outdoor Arthritis Rollator</b>	<b>Antracite grey</b> 312205	<b>Brown metallic</b> 312206
<i>Height of the forearm support (adjustable)</i>	87 - 110 cm	87 - 110 cm
<i>Recommended user heights</i>	140 - 185 cm	140 - 185 cm
<i>Seat height</i>	62 cm	62 cm
<i>Total width opened / folded</i>	71 cm / 31 cm	71 cm / 31 cm
<i>Seat width / Inner width between push handles</i>	56 cm / 46 cm	56 cm / 46 cm
<i>Length (depth)</i>	77 cm	77 cm
<i>Weight (+ bag 380 g)</i>	12 kg	12 kg
<i>Max. load (rollator / bag)</i>	200 kg / 5 kg	200 kg / 5 kg
<i>Brakes</i>	2 lockable with parking function	
<i>Wheel sizes</i>	Front and back: 200 x 35 mm	



# LION ROLLATOR

The Mobilex rollator **Lion** is specially designed for use on off-paved or paved roads. The large soft wheels ensure a comfortable ride even over bumpy cobblestones and uneven gravel or dirt roads. The rollator can easily be folded sideways for a small space requirement and can therefore also be easily transported in the car or parked at home. The rollator comes as standard with a tear resistant, comfortable textile seat, a comfortable and detachable shopping bag, a crutch holder and a back belt. The rollator is available in two metallic colours, brown or anthracite grey, and is tested and approved for a user weight of up to 200 kg.



Lion Rollator	Standard	Small
	312372 (antracite grey) 312374 (metallci brown)	312371 (antracite grey) 312373 (metallic brown)
Height of the push handles (adjustable)	90-112 cm	77-90 cm
Recommended user heights	170 - 220 cm	165 - 175 cm
Seat height	62 cm	53 cm
Total width opened / folded	68 / 26 cm	68 / 26 cm
Inner width between push handles	56 cm	56 cm
Length (depth)	76 cm	76 cm
Weight (+ bag 380 g)	9,3 kg	9,0 kg
Max. load (rollator / bag)	200 kg / 5 kg	200 kg / 5 kg
Brakes	2 lockable with parking function	
Wheel sizes (soft wheels)	Front: 285 x 44 mm - Back: 235 x 44 mm	

# PANTHER ROLLATOR

The Mobilex **Panther** lightweight rollator has a blue aluminium frame. The rollator folds easily sideways and takes up little space when folded, making it easy to transport. The seat is made of durable, long-lasting material that provides comfort and stability for the user. At the front there is a removable shopping bag with two handles. The large, soft wheels make the rollator well suited for both indoor and outdoor use; it is stable and effectively absorbs shocks. The push handles can be adjusted in height very easily and without use of any tools. The rollator is equipped with two lockable brakes that brake on the rear wheels. A tilt pedal is integrated on the rear frame to make it easier to overcome smaller obstacles. The rollator is tested and approved for a user weight of up to 150 kg and is available in two sizes, standard and small.



*A new champagne-colored version of the Panther rollator will be available in 2024 instead of the blue model.*

Panther Rollator	Standard	Small
	3180U1	3180U2
Height of the push handles (adjustable)	78 - 101 cm	71 - 93 cm
Recommended user heights	150 - 195 cm	140 - 187 cm
Seat height	62 cm	55 cm
Total width opened / folded	62 / 22 cm	62 / 22 cm
Inner width between push handles	43 cm	43 cm
Length (depth)	75,5 cm	75,5 cm
Weight (+ bag 380 g)	6,73 kg	6,25 kg
Max. load (rollator / bag)	150 kg / 5 kg	150 kg / 5 kg
Brakes	2 lockable with parking function	
Wheel sizes	Front and back: 200 x 35 mm	



# PUMA AIR ROLLATOR



The **Puma Air** rollator is a high quality outdoor rollator. This rollator has a matte black or red frame and extra large 12" pneumatic tyres (295 mm) for more stability on difficult surfaces. The push handles can easily be adjusted in height without the use of any tools. In addition, the rollator can be folded very easily with a flick of the wrist by pulling the strap on the seat upwards. As standard the Puma Air comes with: a back belt, a cup holder, a crutch holder and a modern, closed shopping bag. With its large pneumatic tyres, the Puma Air is especially suitable for outdoor use, but also for indoor use as it is only 65 cm wide. In total, the outdoor rollator only weighs 9,78 kg and can therefore easily be transported. The rollator is tested and approved for a user weight of up to 150 kg.



<b>Puma Air Rollator</b>	<b>312501+02</b>
Height of the push handles (adjustable)	82 - 99 cm
Recommended user heights	160 - 200 cm
Seat height	59 cm
Total width opened / folded	65 cm / 26 cm
Inner width between push handles	40 cm
Length (depth)	76,4 cm
Weight	9,78 kg
Max. load (rollator / bag)	150 kg / 5 kg
Brakes	2 lockable with parking function
Wheel sizes	Front and back: 295 mm (12") x 45 mm



312560

**Accessory:**  
Tray that is easily attached to the rollator (instead of the back belt).





# JAGUAR ROLLATOR



The **Jaguar** indoor rollator has an aluminium frame and can be easily folded. Overall, the structure is very light, but offers a high level of stability coupled with extreme ease of movement. The rollator is only approved for indoor use and has been tested and approved for a user weight of up to 100 kg.



Rollator Jaguar	Various colours	SUPREME
Height of push handles (adjustable)	83,5 - 94,5 cm	90,3 - 96 cm
Size of the foldet rollator (H x W x D)	83,5 x 54,5 x 27 cm	90,3 x 56,5 x 27 cm
Width	54,5 cm	56,5 cm
Lenght (depth) opened / folded	66 cm / 27 cm	66 cm / 27 cm
Weight	6,3 kg	6,0 kg
Max. load (rollator / bag)	100 kg / 5 kg	100 kg / 5 kg
Brakes	2 lockable with parking function	
Wheel sizes	Front and back: 7" (175 x 27 mm)	

# IMPALA ROLLATOR

The **Impala** rollator has a metallic blue-coloured aluminium frame. The folding system reduces in folded state the space requirements to a minimum. The rollator is equipped with a cane holder, padded seat and backrest and a handy bag under the seat. With its large 200 mm wheels the rollator is suitable for use both outdoors and indoors. The push handles can be easily adjusted in height by the user and provides comfortable and safe walking. The rollator is equipped with two lockable brakes on the rear wheels.




Impala Rollator	Standard 312060	High 312061	Wide 312062
Height of the push handles (adjustable)	79 - 93 cm	90 - 105 cm	79 - 93 cm
Seat height	53 cm	63 cm	53 cm
Total width opened / folded	60 / 27 cm	63 / 35 cm	66 / 27 cm
Inner width between push handles	48 cm	49 cm	53 cm
Length (depth)	70 cm	75 cm	70 cm
Weight	7,4 kg	7,9 kg	7,5 kg
Max. load (rollator / bag)	150 kg / 5 kg	150 kg / 5 kg	150 kg / 5 kg
Brakes	2 lockable with parking function		
Wheel sizes	Front and back: 200 x 35 mm		

# ACCESSORIES

A small selection of accessories available for our Mobilex rollators. The entire selection can be found in our catalog or on [www.mobilex.dk](http://www.mobilex.dk).




Umbrella



One-hand brake



Oxygen holder



Open net bag



Bell



LED lights



High back belt



Tray



Fixed padded seat



Bag with zip



Speed reduction brake

\* Only available for the Buffalo, Gepard, Leopard, Lion and Tiger rollator



# BISON WALKING FRAME

The **Bison** walking frame is a mobility aid for learning and/or training of the active walking and can serve as support in case of balance disturbances. The frame is made of powder-coated aluminium and offers several adjustment functions: the arm rests can be easily adjusted in height by the gas spring support. The handles with the brake levers are angularly adjustable. All these settings can be performed without any use of tools. The manual brakes act on the rear wheels and all four wheels are equipped with parking brakes. The narrow width also allows narrow passages.



<b>Bison walking frame</b>	<b>313300</b>
<i>Width of the frame at the rear wheels</i>	66,0 cm
<i>Width of arm rests min. - max.</i>	66,0 - 80,0 cm
<i>Distance between handles</i>	34,0 cm
<i>Height min. - max.</i>	102,2 - 143,0 cm
<i>Length</i>	86,0 cm
<i>Weight</i>	18,3 kg
<i>Max. load</i>	150 kg

# DYNAWALK WALKING FRAME

The Mobilex **DynaWalk** walking frame is a tool for disabled or elderly people, who need additional support to maintain balance or stability while walking. The frame is made of powder-coated aluminium and offers various adjustment functions: the armrests are adjustable in height, the angle of the handles is adjustable. The safety back bars are easily liftable. The width is easily adjustable through the cross joint. All settings can be done by clips without any use of tools.




<b>DynaWalk walking frame</b>	<b>Standard 313001</b>	<b>Brake 313002</b>
<i>Max. load</i>	135 kg	135 kg
<i>Weight</i>	12,5 kg	12,5 kg
<i>Length</i>	82 cm	82 cm
<i>Width (adjustable)</i>	67,5 - 75 cm	67,5 - 75 cm
<i>Height of armrests (adjustable)</i>	100 - 123 cm	100 - 123 cm
<i>Max. height of the handles</i>	136 cm	136 cm
<i>Front brakes</i>	2 lockable pedal brakes	2 lockable pedal brakes
<i>Rear brakes</i>	2 lockable pedal brakes, with friction	2 lockable brakes on the handles



## MODERN AIDS AT AFFORDABLE PRICES

We offer a wide selection of aids and spare parts for the elderly and people with disabilities. All our products are of the highest quality to meet the special needs of users.



Visit us on [www.mobilex.dk](http://www.mobilex.dk)



Tel. +45 8793 2220  
Fax +45 8793 1777  
[www.mobilex.dk](http://www.mobilex.dk)

Mobilex A/S  
Grønlandsvej 5  
DK-8660 Skanderborg