

KAKADU

COMMODE AND SHOWER CHAIR



Technical data	Small 302015	Standard 302018	Large 302019	XL 302020	Tilt 302024	Tilt L 302025	Diff. with 24"
Seat width	38 cm	45 cm	50 cm	55 cm	45 cm	50 cm	---
Seat depth	38 cm	43 cm	43 cm	43 cm	43 cm	43 cm	---
Seat height	50 - 55 cm	50 - 55 cm	50 - 55 cm	50 - 55 cm	55 - 60 cm	55 - 60 cm	---
Pass over height	42 - 47 cm	42 - 47 cm	42 - 47 cm	42 - 47 cm	---	---	---
Pass over width	29 cm	35 cm	40 cm	45 cm	---	---	---
Total width	49 cm	56 cm	61 cm	66 cm	56 cm	61 cm	+ 8 cm *
Total depth	85 cm	93 cm	93 cm	93 cm	110-120 cm	110-120 cm	+ 5 cm
Total height	96-101 cm	96-101 cm	96-101 cm	96-101 cm	97-141 cm	97-141 cm	---
Max. load	130 Kg	150 Kg	150 Kg	200 kg	130 Kg	130 Kg	---
Weight	11,00 Kg	12,90 Kg	13,50 Kg	14,10 kg	15,10 kg	15,70 kg	+ 2,90 kg
Tilt angle	---	---	---	---	-5° bis +25°	-5° bis +25°	---

* Total width: Difference for Kakadu Tilt (302024 or 302025) with 24" wheels mounted: + 14 cm



The **Kakadu Tilt** commode and shower chairs can be individually adjusted from -5° to +25° in the seat tilt. The headrest with ear recess can be adjusted in height, depth and angle and is included in the delivery. The lower length and the angle of the foot plates are adjustable on the footrests. The complete black upholstery, including the headrest, is made of soft PUR.

The Mobilex **Kakadu** commode and shower chairs are made of high-quality materials and offer the user a high level of comfort. The commode and shower chairs are height adjustable and can be assembled quickly and easily. When developing our Mobilex shower and commode chairs, we have focused on making all models easily adaptable so that they can be used as standard for most users.

