




# RUBBER RAMPS

20 DIFFERENTS HEIGHTS  
IN 3 STANDARD WIDTHS

# RUBBER RAMPS

With rubber threshold ramps, you can easily and quickly equalise the different heights (both inside and outside) that often exist in your home. With rubber ramps installed, you can safely and easily get around your home and access outdoor areas.

- A simple and affordable aid with a long lifespan
- Cost-effective and cost-saving
- With a non-slip surface
- Noise-reducing
- Easy to clean
- For both indoor and outdoor use
- Available with or without glue
- Approved up to a load of 350kg
- Easily customisable and reusable




Available in:

- 20 different heights
- 3 standard widths (700, 900 and 1000 mm)



The ramps are made from recyclable rubber (TPE) and are free of PVC and other harmful phthalates.  
Good for you and the environment.



# ASSEMBLY/MOUNTING

The rubber ramps are available with or without glue on the back:

- On smooth surfaces we recommend rubber ramps *with* glue on the back.
- On sensitive surfaces (e.g. carpets and wooden floors) we recommend rubber ramps *without* glue.


The rubber ramps can easily be shortened in width with a hobby knife to fit your home and your needs.

## Mounting accessories



Item no. 306910

The mounting set for the rubber ramps consists of two mounting brackets that are mounted into the doorstep. The rubber ramps can then be mounted securely and easily without using glue on the surface.



## Adaptation accessories



### Jigsaw blade

Item no. 306904

Suitable for standard jigsaws. Can be used for easy and precise shortening of all Mobilex rubber ramps.



### Double adhesive tape

Item no. 306911 / 306912

If you want to use your ramp elsewhere, you can use this adhesive tape. Available in two widths: 20 and 35 mm.



### Degreaser for glue removal


Item no. 306905

Remove glue remains on either rubber ramps or surfaces with this cleaner. The bottle contains 0.5 L.

## WE ALSO HAVE A WIDE RANGE OF TRANSPORTABLE ALUMINIUM WHEELCHAIR RAMPS:



SC ramps



Foldable doorstep ramps



Foldable CH channel ramps

## RUBBER RAMPS WITH GLUE

RUBBER RAMPS 700 MM	
Item no.	Size (H x D x W)
30610L	10 x 100 x 700 mm
30612L	12 x 120 x 700 mm
30616L	16 x 150 x 700 mm
30620L	20 x 150 x 700 mm
30624L	24 x 150 x 700 mm
30628L	28 x 200 x 700 mm
30632L	32 x 200 x 700 mm
30636L	36 x 200 x 700 mm
30640L	40 x 250 x 700 mm
30644L	44 x 250 x 700 mm
30648L	48 x 250 x 700 mm

RUBBER RAMPS 900 MM	
Item no.	Size (H x D x W)
306040	4 x 40 x 900 mm
306060	6 x 60 x 900 mm
306080	8 x 80 x 900 mm
306100	10 x 100 x 900 mm
306120	12 x 120 x 900 mm
306160	16 x 150 x 900 mm
306200	20 x 150 x 900 mm
306240	24 x 150 x 900 mm
306280	28 x 200 x 900 mm
306320	32 x 200 x 900 mm
306360	36 x 200 x 900 mm
306400	40 x 250 x 900 mm
306440	44 x 250 x 900 mm
306480	48 x 250 x 900 mm
306520	52 x 300 x 900 mm
306560	56 x 300 x 900 mm
306600	60 x 300 x 900 mm
306640	64 x 350 x 900 mm
306680	68 x 350 x 900 mm
306720	72 x 350 x 900 mm

RUBBER RAMPS 1000 MM	
Item no.	Size (H x D x W)
306042	4 x 40 x 1000 mm
306062	6 x 60 x 1000 mm
306082	8 x 80 x 1000 mm
306102	10 x 100 x 1000 mm
306122	12 x 120 x 1000 mm
306162	16 x 150 x 1000 mm
306202	20 x 150 x 1000 mm
306242	24 x 150 x 1000 mm
306283	28 x 200 x 1000 mm

## RUBBER RAMPS WITHOUT GLUE

RUBBER RAMPS 700 MM	
Item no.	Size (H x D x W)
30652M	52 x 300 x 700 mm
30656M	56 x 300 x 700 mm
30660M	60 x 300 x 700 mm
30664M	64 x 350 x 700 mm
30668M	68 x 350 x 700 mm
30672M	72 x 350 x 700 mm

RUBBER RAMPS 900 MM	
Item no.	Size (H x D x W)
306041	4 x 40 x 900 mm
306061	6 x 60 x 900 mm
306081	8 x 80 x 900 mm
306101	10 x 100 x 900 mm
306121	12 x 120 x 900 mm
306161	16 x 150 x 900 mm
306201	20 x 150 x 900 mm
306241	24 x 150 x 900 mm
306281	28 x 200 x 900 mm
306321	32 x 200 x 900 mm
306361	36 x 200 x 900 mm
306401	40 x 250 x 900 mm
306441	44 x 250 x 900 mm
306481	48 x 250 x 900 mm
306521	52 x 300 x 900 mm
306561	56 x 300 x 900 mm
306601	60 x 300 x 900 mm
306641	64 x 350 x 900 mm
306681	68 x 350 x 900 mm
306721	72 x 350 x 900 mm

RUBBER RAMPS 1000 MM	
Item no.	Size (H x D x W)
306282	28 x 200 x 1000 mm
306322	32 x 200 x 1000 mm
306362	36 x 200 x 1000 mm
306402	40 x 250 x 1000 mm
306442	44 x 250 x 1000 mm
306482	48 x 250 x 1000 mm
306522	52 x 300 x 1000 mm
306562	56 x 300 x 1000 mm
306602	60 x 300 x 1000 mm
306642	64 x 350 x 1000 mm
306682	68 x 350 x 1000 mm
306722	72 x 350 x 1000 mm

Find and order all rubber ramps directly at [www.mobilex.dk](http://www.mobilex.dk)

